LISA - Lumbar Implant for Stiffness Augmentation

EN
Instructions For Use - Implants
LISA Posterior Stabilization System

IMPORTANT INFORMATION ABOUT LISA - Lumbar Implant for Stiffness Augmentation (IMPLANT)
 Before using a product sold by BackBone, surgeons should carefully study the recommendations, warnings and instructions, as well as specific information available for the particular product line (Instruments for use, surgical technique, etc...). BackBone shall not be liable for any complications resulting from the use of this product line outside of its specified areas of use, its indications, contraindications and recommendations given in the Instructions For Use and/or the surgical technique.

Compatibility has not been established between all of the product lines. Only authorized combinations of products should be used. Devices must be used only as indicated in the Instructions For Use and/or the surgical technique.

1. DESCRIPTION
 The LISA implant is intended to stabilize the spinal column by using an operative technique referred to as "non-fusion". The LISA implant is inserted by posterior approach and dedicated for the lumbar area (from L1 to L5).

The implant is comprised of a spacer designed to be implanted between two lumbar spinous processes and available in various sizes. This spacer is provided with a polyester braud used to hook the spinous processes and a hooker made of titanium alloy to lock the braud into the spacer. The LISA implant shape is adapted to the morphology of the intervertebral spaces of the treated levels, as well as to the recommended surgical methods.

For implanting the system, please refer to the **LISA - Surgical Technique**.

2. INDICATIONS
 The LISA Posterior Dynamic Stabilization System treats low-back pain that accompanies degenerative lesions of grade I, III and IV (Pfirrmann MRI classification).

3. CONTRAINDICATIONS
 a. Stage V degenerative disk lesions in Pfirrmann's MRI classification.
 b. Spondylosisthesis.
 c. Osteoporosis.
 d. Non-specific back pain.
 e. Modic 2 and Modic 3 changes.
 f. This implant is not indicated for the L5/S1 segments.
 g. Local or general infections that may compromise the surgical goals.
 h. Major local inflammatory phenomena.

i. Pregnancy.
 j. Immunosuppressive diseases.
 k. Bone immaturity.
 l. Severe mental illnesses.
 m. Bone metabolism diseases that may compromise the mechanical support expected from this type of implant.
 n. Excessive physical activities.

4. POSSIBLE SIDE EFFECTS
 All potential side effects associated with spinal surgery independent of the medical device are possible. With the use of this medical device, the list of potential side effects includes, among others:

- Neurological complications, paralysis, soft tissue injuries, pain, implant migration, erosion or implant breakage.
- Infection or deep infections and inflammatory phenomena.
- Allergic reactions to the materials comprising the implant.
- Abration of the bone density due to a change in the distribution of mechanical stresses.
- Neurological injuries and/or injuries to the dura mater during the surgical procedure.
- Though the risks have been justified at a low level, the use of Magnetic Resonance Imaging (MRI) on patients treated with LISA implants may result in possible side effects such as migration or localized heat generation due to the metallic component of the LISA implant (Locker).

5. WARNINGS
 The LISA implant is a single use device and the reuse of LISA may cause infections or ineffective cures.

These devices must be never re-sterilized. Potential risks related to re-sterilization of the device that might affect the patient health and safety include:

- The transmission of infectious or viral agents. The implant may not be re-cleaned or re-sterilized.
- Loss of the functional and mechanical properties of the implant (including possible rupture) after the first implantation and subsequent removal of the device.

Patient physical activity should be handled with caution in order to avoid complications of the device that might affect the patient health and safety include:

- The transmission of infectious or viral agents: no re-sterilization method has been validated for this device.
- Change in the physical properties of the material composing the device that lead to loss of its functionality and/or injuries to the dura mater during the surgical procedure.
- Any contaminated implant should be treated as biological waste.
- The implant may impede localized medical procedures such as lumbar punctures or spinal anesthesia.

6. MATERIALS
 LISA implants constituent materials are:
 - PEEK (poly-ether-ether-ketone) polymer
 - Ti6Al4V titanium alloy (ISO 5832/3)
 - Woven PET (polyethylene-terephthalate)

7. PRE-OPERATIVE PRECAUTIONS

a. Patient's weight: overweight conditions cause additional stresses that may lead, in combination with other factors, to rupture of the implants.
 b. Mental handicap: there are a greater risks in patients who are unable to follow the surgeon's recommendations.

c. Hypersensitivity to PEEK and/or PET and/or consistent metals: if hypersensitivity is suspected or confirmed, it is recommended that the patient's tolerance of the substances comprising the material be checked prior to inserting the device.

8. PER-OPERATIVE PRECAUTIONS
 The details of the operative instructions are found in the LISA Surgical Technique supplied by BackBone.

a. Insertion of an implant must be done using the instruments designed and supplied for this purpose and the specific technique to each device.
 b. Bone quality: a case of osteoporosis, or any other tissue disease that may alter the mechanical properties of the spinous processes must be taken into consideration when deciding to use a LISA implant.
 c. It is imperative that the level of tension given by the tensioner used simultaneously with the torque limiting handpiece be followed. If the user overtightens beyond the recommended tension, the resulting tension on the braud may damage the spinous processes, depending on the patient's bone quality.

9. POST-OPERATIVE PRECAUTIONS

a. The surgeon should warn the patient about the precautions to be taken after implantation of the device.
 b. A rigid external lumbar support is normally not required. However, this decision is up to the surgeon, depending on each patient (bone quality, treated and related diseases, patient level of activity and weight, etc...) .
 c. Patient physical activity: intense physical activity increases the risk of mobility, deformation and rupture of the implants.
 d. A physical handicap will require special attention or adaptation of the post-operative rehabilitation method.

10. PACKAGING
 The implants are delivered in sterile packaging, which must be intact upon receipt and before its use.

11. RECOMMENDATIONS FOR IMPLANTS SUPPLIED IN STERILE CONDITION
 It is recommended that the operator should be given clear advice. BackBone shall not be held liable for any use of products after their expiration date. Sterile implants should not be allowed to come into contact with disinfectant solutions or products.
 It is recommended that the integrity of the packaging be verified before use. Sterility is only valid if the packaging shows no signs of deterioration. If the packaging is damaged, or after it is opened, the implant must not be re-sterilized under any circumstance.

12. ADDITIONAL INFORMATION
 The LISA – Surgical Technique is available upon request from your BackBone representative or directly from BackBone. If the Surgical Technique in the user's possession is older than two years prior to the date of surgery, it is recommended that you request an updated version.

13. COMPLAINTS
 Any healthcare professional (customer or user) who has a complaint concerning the services and/or quality, identification, robustness, reliability, safety, efficacy and/or performance of BackBone products should report it to their representative or distributor of these products. The distributor will inform BackBone of this complaint as soon as possible, using an incident report form. The company should be immediately alerted by phone or by any written means in the case of a malfunction or deterioration of the device, or of any inappropriate part of the Instructions for Use that has caused, or might have caused, death or serious deterioration in the health of a patient or user.

Medical device vigilance: incident report forms must contain as many details as possible, such as identification of the product (description, batch number, reference number), the type of complaint or a detailed description of the incident, as well as its consequences. Incident report forms must be submitted along with any technical aspects that might facilitate an expert appraisal, such as the implant itself, X-rays, etc. . . .

12. SUPPLEMENTARY INFORMATION
 The LISA – Surgical Technique (Kinurg) can be found at www.backbone.fr or you can contact BackBone representative directly via BackBone. This new update brings us a further technical update of the kinurg for both L1 and L5.

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