## Instructions for use LISA Implants

BB-LISA-8-001-V03 March 2019

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### EN Instructions for Use - Implants LISA Posterior Stabilization System

#### IMPORTANT INFORMATION ABOUT LISA - Lumbar Implant for Stiffness Augmentation (IMPLANT)

Before using a product sold by BackBone, surgeons should carefully study the recommendations, warnings and instructions, as well as specific information available for the particular product line (instructions for use, surgical technique, etc. ...). BackBone shall not be liable for any complications resulting from the use of this product line outside of its specified areas of use, its indications, contraindications and recommendations given in the Instructions For Use and/or the surgical technique.

Compatibility has not been established between all of the product lines. Only authorized combinations of products should be used. Devices must be used only as indicated in the Instructions For Use and/or the surgical technique.

**1. DESCRIPTION**
The LISA implant is intended to stabilize the spinal column by using an operative technique referred to as "non-fusion". The LISA Implant is inserted by posterior approach and dedicated for the lumbar area (from L1 to L5).
The implant is composed of a spacer designed to be implanted between two lumbar spinous processes and available in various sizes. This spacer is provided with a polyester braid used to hook the spinous processes and a blocker made in titanium alloy to lock the braid into the spacer.
The LISA Implant shape is adapted to the morphology of the intervertebral spaces of the treated levels, as well as to the recommended surgical methods.
For implanting the system, please refer to the LISA - Surgical Technique.

**2. INDICATIONS**
The LISA Posterior Dynamic Stabilization System treats low-back pain that accompanies degenerative lesions of grade II, III and IV (Pfirrmann MRI classification).

- 3. CONTRAINDICATIONS**
- Stage V degenerative disk lesions in Pfirmann's MRI classification.
  - Spondylolithesis.
  - Osteoporosis.
  - Non-specific back pain.
  - Modic 2 and Modic 3 changes.
  - This implant is not indicated for the L5/S1 segments.
  - Local or general infections that may compromise the surgical goals.
  - Major local inflammatory phenomena.
  - Pregnancy.
  - Immunosuppressive diseases.
  - Bone immaturity.
  - Severe mental illnesses.
  - Bone metabolism diseases that may compromise the mechanical support expected from this type of implant.
  - Excessive physical activities.

#### 4. POSSIBLE SIDE EFFECTS

All potential side effects associated with spinal surgery independent of the medical device are possible. With the use of this medical device, the list of potential side effects includes, among other:

- Neurological complications, paralysis, soft tissue injuries, pain, implant migration, erosion or implant breakage.
- Neurologic or deep infections and inflammatory phenomena.
- Allergic reactions to the materials composing the implant.
- Alteration of the bone density due to a change in the distribution of mechanical stresses.
- Neurological injuries and/or injuries to the dura mater during the surgical procedure.
- Though the risks have been judged at a low level, the use of Magnetic Resonance Imaging (MRI) on patients treated with LISA implants may result in possible side effects such as migration or localized heat generation due to the metallic component of the LISA implant (Locker).

#### 5. WARNINGS

- The LISA implant is a single use device and the reuse of LISA may cause infections or ineffective cure.
- These devices must be implanted by surgeons who have been properly trained in spinal surgery. The decision to implant them should be made only after taking into consideration the medical and surgical indications, contraindications, side effects and precautions contained in these Instructions For Use and the limitations of this type of surgery.
- Even if a device seems to be intact after it has been removed from a patient, these implants should never be re-used. Potential risks related to re-use of the device that might affect a patient's health and safety include:
  - The transmission of infectious or viral agents. The implant may not be re-cleaned or re-sterilized.
  - Loss of the functional and mechanical properties of the implant (including possible rupture) after the first implantation and subsequent removal of the device.
- Sterile Implants must be never re-sterilized. Potential risks related to re-sterilization of the device that might affect the patient health and safety include:
  - The transmission of infectious or viral agents; no re-sterilization method has been validated for this device.
  - Change in the physical properties of the material composing the device that lead to loss of functionality and mechanical properties, including rupture or degradation of the device.
- Any contaminated implant should be treated as biological waste
- The implant may impede localized medical procedures such as lumbar punctures or spinal anesthesia.

#### 6. MATERIALS

- LISA Implants constituent materials are:
- PEEK (poly-ether-ether-ketone) polymer
  - Ti6Al4V titanium alloy (ISO 5832/3)
  - Woven PET (polyethylene terephthalate)

#### 7. PRE-OPERATIVE PRECAUTIONS

- Patient's weight: overweight conditions cause additional stresses that may lead, in combination with other factors, to rupture of the implants.
- Mental handicap: there are a greater risks in patients who are unable to follow the surgeon's recommendations.
- Hypersensitivity to PEEK and/or PET and/or constituent metals: If hypersensitivity is suspected or confirmed, it is recommended that the patient's tolerance of the substances comprising the implant be checked prior to inserting the device.

#### 8. PER-OPERATIVE PRECAUTIONS

The details of the operative instructions are found in the LISA Surgical Technique supplied by BackBone.

- Insertion of an implant must be done using the instruments designed and supplied for this purpose and the specific technique to each device.
- Bone quality: a case of osteoporosis, or any other tissue diseases that may alter the mechanical properties of the spinous processes must be taken into consideration when deciding to use a LISA implant.
- It is imperative that the level of tension given by the tensioner used simultaneously with the torque limiting handle be followed. If the user overtightens beyond the recommended tension, the resulting tension on the band may damage the spinous processes, depending on the patient's bone quality.

#### 9. POST-OPERATIVE PRECAUTIONS

- The surgeon should warn the patient about the precautions to be taken after implantation of the device.
- A rigid external lumbar support is normally not required. However, this decision is up to the surgeon, depending on each patient (bone quality, treated and related diseases, patient level of activity and weight, etc. ...).
- Patient physical activity: intense physical activity increases the risk of mobility, deformation and rupture of the implants.
- A physical handicap will require special attention or adaptation of the post-operative rehabilitation method.

#### 10. PACKAGING

The implants are delivered in sterile packaging, which must be intact upon receipt and before its use.

#### 11. RECOMMENDATIONS FOR IMPLANTS SUPPLIED IN STERILE CONDITION

- It is recommended that the expiration date be checked prior to use. BackBone shall not be held liable for any use of products after their expiration date. Sterile implants should not be allowed to come into contact with disinfectant solutions or products.
- It is recommended that the integrity of the packaging be verified before use. Sterility is only valid if the packaging shows no signs of deterioration. If the packaging is damaged, or after it is opened, the implant must not be re-sterilized under any circumstance.

#### 12. ADDITIONAL INFORMATION

The LISA – Surgical Technique is available upon request from your BackBone representative or directly from BackBone. If the Surgical Technique in your s' possession is older than two years prior to the date of surgery, it is recommended that you request an updated version.

#### 13. COMPLAINTS

Any healthcare professional (customer or user) who has a complaint concerning the services and/or quality, identification, robustness, reliability, safety, efficacy and/or performance of BackBone products should report it to their representative or distributor of these products. The distributor will inform BackBone of this complaint as soon as possible, using an incident report form. The company should be immediately alerted by phone or by any written means in the case of a malfunction or deterioration of the device, or of any inappropriate part of the Instructions For Use that has caused, or might have caused, death or serious deterioration in the health of a patient or user.

Medical device vigilance: Incident report forms must contain as many details as possible, such as identification of the product (description, batch number, reference number), the type of complaint or a detailed description of the incident, as well as its consequences. Incident report forms must be submitted along with any technical aspects that might facilitate an expert appraisal, such as the implant itself, X-rays, etc. ...

For further information or to submit a complaint, contact:

**BACKBONE**  
81, Boulevard Pierre Ter  
33110 LE BOUSCAT  
France  
Tel. +33 (0) 556 424 430

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