Instructions For Use - Implants LISA Posterior Stabilization System

IMPORTANT INFORMATION ABOUT LISA - Lumbar Implant for Stiffness Augmentation (IMPLANT)

Before using a product sold by BACKBONE, surgeons should carefully study the recommendations, warnings, and instructions, as well as specific information available for the particular product line (instructions for use, surgical technique, etc.) BACKBONE shall not be liable for any complications resulting from the use of this product line outside of its specified areas of use, its indications, contraindications, and recommendations given in the Instructions For Use, and/or the surgical technique.

Compatibility has been established between LISA Implants and LISA Instruments. Only authorized combinations of products between LISA System should be used. Devices must be used only as indicated in the Instructions For Use and/or the surgical technique.

The implant is placed by the associated LISA instruments

For implanting the system, please refer to the **LISA - Surgical Technique**.

LISA device is a posterior lumbar dynamic stabilization system designed to stabilize the treated level while preserving motion.

The LISA device consists of 3 components packaged and sterilized separately: A PEEK interspinous spacer, a polyester band, and a titanium blocker. The spacer is positioned between two adjacent spinous processes, the band is belted around the spinous processes and through the spacer, and the blocker is used to lock the band inside the spacer.

1. INTENDED PURPOSE

The intended purpose of the LISA device is to safely improve back pain, leg pain and disability while allowing motion preservation between two adjacent lumbar vertebrae when used in degenerative lesions of grade II, III, IV according to Pfirrmann MRI classification. It can be used in up to two adjacent levels from L1 to L5.

1.1. INTENDED USERS

LISA devices must be implanted by surgeons who have been properly trained in spinal surgery. The decision to implant them should be made only after taking into consideration the medical and surgical indications, contraindications, side effects and precautions contained in the Instructions For Use and the limitations of this type of surgery.

1.2. INTENDED TARGET POPULATION

LISA is intended to be used on skeletally mature patients suffering from low-back pain that accompanies degenerative lesions of grade II, III and IV (Pfirrmann MRI classification), in accordance with the indications and contra-indications of the device.

1.3. INDICATIONS

The LISA Posterior Dynamic Stabilization System treats low-back pain that accompanies degenerative lesions of grades II, III, and IV (Pfirrmann MRI classification).

1.4. CONTRAINDICATIONS

Contraindications:

- a. Stage V degenerative disk lesions in Pfirrmann's MRI classification.
- b. Spondylolisthesis.

- c. Osteoporosis.
- d. Non-specific back pain.
- e. Modic 2 and Modic 3 changes.
- f. This device is not indicated for the L5/S1 segments.
- g. Local or general infections that may compromise the surgical goals.
- h. Major local inflammatory phenomena.
- i. Pregnancy.
- j. Immunosuppressive diseases.
- k. Bone immaturity.
- l. Severe mental illnesses.
- m. Bone metabolism diseases that may compromise the mechanical support expected from this type of implant.
- n. Excessive physical activities.

1.5. WARNINGS

- The LISA implant is a single-use device, and the reuse of LISA may cause infections or ineffective care.
- Sterile implants must never be re-sterilized. Potential risks related to re-sterilization of the device that might affect the patient health and safety include:
 - The transmission of infectious or viral agents: no re-sterilization method has been validated for this device.
 - Change in the physical properties of the material composing the device leads to loss of functionality and mechanical properties, including rupture or degradation of the device.
- Even if a device seems intact after being removed from a patient, these implants should never be re-used. Potential risks related to the re-use of the device that might affect a patient's health and safety include:
 - The transmission of infectious or viral agents. The implant may not be re-cleaned or re-sterilized.
 - Loss of the functional and mechanical properties of the implant (including possible rupture) after the first implantation and subsequent removal of the device.
- Any contaminated implant should be treated as biological waste.
- The implant may impede localized medical procedures such as lumbar punctures or spinal anesthesia.

1.6. PRECAUTIONS

1.6.1. PRE-OPERATIVE PRECAUTIONS

- a. Patient's weight: overweight conditions cause additional stresses that may lead, in combination with other factors, to the rupture of the implants.
- b. Mental handicap: there is a greater risk in patients who cannot follow the surgeon's recommendations.
- c. Hypersensitivity to PEEK and/or PET and/or constituent metals: if hypersensitivity is suspected or confirmed, it is recommended that the patient's tolerance of the substances comprising the implant be checked before inserting the device.

1.6.2. PER-OPERATIVE PRECAUTIONS

The details of the operative instructions are found in the LISA Surgical Technique supplied by BACKBONE.

- a. Insertion of an implant must be done using the instruments designed and supplied for this purpose and the specific technique for each device.
- b. Bone quality: a case of osteoporosis or any other tissue disease that may alter the spinous processes' mechanical properties must be considered when deciding to use a **LISA** implant.

- c. It is imperative that the level of tension given by the tensioner used simultaneously with the torque limiting handle be followed. If the user overtightens beyond the recommended tension, the resulting tension on the band may damage the spinous processes, depending on the patient's bone quality.

1.6.3. POST-OPERATIVE PRECAUTIONS

The surgeon should warn the patient about the precautions to be taken after the implantation of the device. If the performance of the device changes from what the surgeon indicated, the patient must contact the surgeon.

- a. A rigid external lumbar support usually is not required. However, this decision is up to the surgeon, depending on each patient (bone quality, treated and related diseases, patient-level of activity and weight, etc...).
- b. Patient physical activity: intense physical activity increases the risk of mobility, deformation and rupture of the implants.
- c. A physical handicap will require special attention or adaptation to the post-operative rehabilitation method.
- d. After the implantation of LISA, the surgeon gives to the patient the implant card completed with the identification labels of the LISA implants used.

1.7. ADVERSE EFFECTS

All potential adverse effects of spinal surgery independent of the medical device are possible. The adverse effects include, among others:

- Neurological complications, paralysis, soft tissue injuries, pain,
- Superficial or deep infections and inflammatory phenomena
- Spinous Process Fractures
- Herniated disc/Recurrence of herniated disc
- Residual stenosis
- Neurological injuries and/or damages to the dura mater during the surgical procedure
- Alteration of the bone density due to a change in the distribution of mechanical stresses

With the use of implants from the LISA dynamic stabilization system, the list of potential adverse effect may include:

- Device migration, dislodgment, implant loosening or breakage.
- Spinous Process Fractures
- Allergic reactions to the materials comprising the implant.
- Heating or migration of the implant following the use of magnetic resonance imaging
- Neurological complications following the device use
- Paralysis following the device use following the device use
- Though the pain is reduced, the pain is not sufficiently contained following LISA implantation
- Superficial or deep infections following the device use
- Inflammatory phenomena following the device use
- Alteration of the bone density due to a change in the distribution of mechanical stresses following the device use
- Duramater injury following the device use
- New stenosis after the use of LISA
- Adjacent level slip
- Modic changes in endplate due to the LISA implantation
- Recurrent disc herniation due to the LISA implantation

1.8. RESIDUAL RISKS

1. Torque Limiting Handle: The torque limiting handle is required to be used to limit the tightening of the band around the spinous processes. With an excessive tightening, there is a risk of spinous process fracture during the surgery of short term after the surgery. The limiting handle and the

torque limitation have been defined based on the literature, and the device has been designed and produced to verify that the torque limiting handle achieves its performance. The handle is required to avoid the risk of spinous process fracture.

2. Leachable substance from the band inside the patient: The band raw materials have been selected to be compatible with patient safety. All tests complied with the acceptance criteria and met the expectations of the respective standards though a slight irritation was observed. Possible side effects may include allergic reactions to materials of the implant and inflammatory phenomenon.

2. CLINICAL BENEFITS

The following clinical benefits are expected for patients operated with LISA:

- Reduction of disability in daily activities after operation (post-operative vs. pre-operative) (e.g. pain intensity, personal care, lifting, walking, etc.)
- Back pain reduction (post-operative vs. pre-operative)
- Leg pain reduction (post-operative vs. pre-operative)
- Blood loss during LISA surgery acceptable when compared to blood loss during a surgery with decompression and fusion (alternative to LISA)
- Surgery length for LISA surgery acceptable when compared to surgery length for a surgery with decompression and fusion (alternative to LISA)
- Number of days at hospital after LISA operation acceptable when compared to number of days at hospital after a surgery with decompression and fusion (alternative to LISA)

3. PERFORMANCE CHARACTERISTICS OF THE DEVICE

The following performance characteristics are expected for patients operated with LISA:

- Mobility preservation after surgery (in terms of Range of Motion) at the operated and adjacent levels.
- Protection of adjacent levels from degeneration

4. MATERIALS

LISA implants constituent materials are:

- The spacer is made of PEEK (poly-ether-ether-ketone) polymer. It is in contact with spinous process, blood and soft tissue. The volume of substance-exposed to the patient is :
 - o Spacer size 6 : $V = 3153,9 \text{ mm}^3$
 - o Spacer size 8 : $V = 3715,38 \text{ mm}^3$
 - o Spacer size 10 : $V = 4337,71 \text{ mm}^3$
 - o Spacer size 12 : $V = 4917,97 \text{ mm}^3$
- The blocker is made of Ti6Al4V titanium alloy (ISO 5832/3). It is in contact with blood and soft tissue. The volume of substance-exposed to the patient is $V = 312.24 \text{ mm}^3$.
- The band is made of Woven PET (polyethylene terephthalate). It is in contact with the spinous process, blood, and soft tissue. The surface of substance-exposed to the patient is $S = 10253 \text{ mm}^2$.

Characteristics of the devices are resumed in the Table below:

LISA	Material	Surface implant (mm ²)	Implant volume (mm ³)	Mass (gr)
Band	PET (polyethylene terephthalate)	10253	-	3.42

Blocker	Ti6Al4V titanium alloy (ISO 5832-3)	423	312	1,38
Size Spacer 6	PEEK (poly-ether-ether-ketone) polymer	2488	3154	4,16
Size Spacer 8	PEEK (poly-ether-ether-ketone) polymer	2732	3715	4,90
Size Spacer 10	PEEK (poly-ether-ether-ketone) polymer	2975	4338	5,72
Size Spacer 12	PEEK (poly-ether-ether-ketone) polymer	3214	4918	6,49

5. PACKAGING

The implants are delivered in sterile packaging, which must be intact upon receipt and before use.

The following information is available on the LISA Implants label:

- Medical device
- The trade name in different languages
- Reference of the product
- Address of the manufacturer
- Batch number
- UDI-DI
- UDI Data Matrix
- Expiration date
- CE marking and BSI number
- e-IFU information and BACKBONE website
- LISA Implant cannot be re-sterilized
- LISA Implant is a single-use medical device
- Sterility indicator (must be red if exposed to gamma sterilization)
- MR conditional
- Single sterile barrier system with protective packaging inside

6. RECOMMENDATIONS FOR IMPLANTS SUPPLIED IN STERILE CONDITION

The device is supplied in a sterile packaging gamma irradiated.

It is recommended that the expiration date be checked before use. BACKBONE shall not be held liable for any use of products after their expiration date. Sterile implants should not be allowed to come into contact with disinfectant solutions or products.

It is recommended that the integrity of the packaging be verified before use. Sterility is only valid if the packaging shows no signs of deterioration. If the packaging is damaged or after it is opened, the implant must not be implanted, not be re-sterilized under any circumstance.

The sterility indicator placed on the label must be red before use. The product cannot be used in if the sterility indicator is green (not sterilized).

No special storage and environmental conditions are required for the LISA implant. The LISA implant must be handled by hospital staff and surgeons who have been appropriately trained with the LISA implant and the instruments.

7. MRI SAFETY INFORMATION

A patient with LISA implant may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

Parameter	Condition of Use / Information
Static Magnetic Field Strength (B0)	1,5 T and 3 T

BO Field Orientation	Horizontal, Cylindrical Bore
Maximum Spatial Field Gradient	1.5 T: 82.21 T/m (8221 G/cm) 3 T: 41.11 T/m (4111 G/cm)
RF Polarization (Note: formerly called RF Excitation)	1.5 T: Circularly Polarized (CP) 3 T: Multichannel-2 (MC-2) or Circularly Polarized (CP)
RF Transmit Coil	Whole Body transmit RF coil or extremities local transmit RF coils
RF Receive Coil	Any receive RF coil may be used
MR System (RF) Operating Modes or Constraints	Normal Operating Mode
Whole Body Averaged SAR	Whole Body Averaged SAR \leq 2 W/kg
Item Configuration	Up to 2 LISA devices implanted in a patient
Scan Duration and Wait Time	Maximum 1 hour of continuous scanning without a cooling period
MR Image Artifact	The presence of this implant within the MR image field of view may produce an image artifact. Some manipulation of scan parameters may be needed to compensate for the artifact. In non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the device edge when imaged with a spin echo pulse sequence and 14 mm with a gradient echo, both at 1.5 T.

8. ADDITIONAL INFORMATION:

- The LISA – Surgical Technique is available upon request from your BACKBONE representative or directly from BACKBONE. If the Surgical Technique in the user’s possession is older than two years before the date of surgery, it is recommended that you request an updated version.

A summary of the clinical data and information regarding the device's safety and clinical performance (SSCP) is available on the website: <https://backbone.pro/patients-information/>.

The SSCP is also available on EUDAMED. To find it, please use one of the following Basic UDI-DI:

- LISA Band: 376024863LISA101FT
- LISA Blocker: 376024863LISA104FZ
- LISA Spacer: 376024863LISA106G5

9. COMPLAINTS

Any healthcare professional (customer or user) who has a complaint concerning the services and/or quality, identification, robustness, reliability, safety, efficacy, and/or performance of BACKBONE products should report it to their representative or distributor of these products. Using an incident report form, the distributor will inform BACKBONE of this complaint as soon as possible. The company should be immediately alerted by phone or by any written means in the case of a malfunction or deterioration of the device (eg. packaging) or any inappropriate part of the Instructions For Use that has caused or might have caused death or serious deterioration in the health of a patient or user.

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Medical device vigilance: Incident report forms must contain as many details as possible, such as identification of the product (description, batch number, reference number), the type of complaint, and a detailed description of the incident, as well as its consequences. Incident report forms must be submitted along with any technical aspects that might facilitate an expert appraisal, such as the device

itself, X-rays, etc.

For further information or to submit a complaint, please contact:

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