	SSCP LISA - Lumbar Implant for Stiffness Augmentation	BF-127-MNGQ-V01
		Datum verzije: 24/01/2022

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9. Bibliografske reference

1. Nicholson JA, Scott CEH, Duckworth AD, Burke JG, Gibson JNA. Survival analysis of the Wallis interspinous spacer used as an augment to lumbar decompression. *Br J Neurosurg.* déc 2017;31(6):688-94.
2. Boody BS, Smucker JD, Sasso W, Miller JW, Snowden R, Sasso RC. Evaluation of DIAM™ Spinal Stabilization System for lower lumbar disc degenerative disease: A randomized, prospective, single-site study. *J Orthop.* oct 2020;21:171-7.
3. Mo Z, Li D, Zhang R, Chang M, Yang B, Tang S. Comparative effectiveness and safety of posterior lumbar interbody fusion, Coflex, Wallis, and X-stop for lumbar degenerative diseases: A systematic review and network meta-analysis. *Clin Neurol Neurosurg.* sept 2018;172:74-81.
4. Deer TR, Grider JS, Pope JE, Falowski S, Lamer TJ, Calodney A, et al. The MIST Guidelines: The Lumbar Spinal Stenosis Consensus Group Guidelines for Minimally Invasive Spine Treatment. *Pain Pract Off J World Inst Pain.* mars 2019;19(3):250-74.
5. Zhou Z, Jin X, Wang C, Wang L. Wallis interspinous device versus discectomy for lumbar disc herniation : A comparative study. *Orthopade.* févr 2019;48(2):165-9.
6. Pintauro M, Duffy A, Vahedi P, Rymarczuk G, Heller J. Interspinous implants: are the new implants better than the last generation? A review. *Curr Rev Musculoskelet Med.* juin 2017;10(2):189-98.
7. Pham MH, Mehta VA, Patel NN, Jakoi AM, Hsieh PC, Liu JC, et al. Complications associated with the Dynesys dynamic stabilization system: a comprehensive review of the literature. *Neurosurg Focus.* janv 2016;40(1):E2.
8. Cho HJ, Ko YS, Won YI, Lee CH, Yang SH, Kim CH, et al. The Efficacy of Lumbar Hybrid Fusion for the Prevention of Adjacent Segment Disease: Fact or Artifact? A Meta-analysis. *Clin Spine Surg.* 9 oct 2020;
9. Porrino J, Rao A, Moran J, Wang A, Grauer J, Haims A, et al. Current concepts of spondylosis and posterior spinal motion preservation for radiologists. *Skeletal Radiol.* nov 2021;50(11):2169-84.
10. Chopko B, Liu JC, Khan MK. Anatomic surgical management of chronic low back pain. *Neuromodulation J Int Neuromodulation Soc.* oct 2014;17 Suppl 2:46-51.
11. Wang K, Zhu Z, Wang B, Zhu Y, Liu H. Bone resorption during the first year after implantation of a single-segment dynamic interspinous stabilization device and its risk factors. *BMC Musculoskelet Disord.* 14 mai 2015;16(1):117.
12. Krappel F, Brayda-Bruno M, Alessi G, Remacle JM, Lopez LA, Fernández JJ, et al. Herniectomy versus herniectomy with the DIAM spinal stabilization system in patients with sciatica and concomitant low back pain: results of a prospective randomized controlled multicenter trial. *Eur Spine J Off Publ Eur Spine Soc Eur Spinal Deform Soc Eur Sect Cerv Spine Res Soc.* mars 2017;26(3):865-76.
13. Korovessis P, Repantis T, Zacharatos S, Zafiroopoulos A. Does Wallis implant reduce adjacent segment degeneration above lumbosacral instrumented fusion? *Eur Spine J Off Publ Eur Spine Soc Eur Spinal Deform Soc Eur Sect Cerv Spine Res Soc.* juin 2009;18(6):830-40.
14. Lee SH, Seol A, Cho TY, Kim SY, Kim DJ, Lim HM. A Systematic Review of Interspinous Dynamic Stabilization. *Clin Orthop Surg.* sept 2015;7(3):323-9.
15. Sénégas J. Mechanical supplementation by non-rigid fixation in degenerative intervertebral lumbar segments: the Wallis system. *Eur Spine J Off Publ Eur Spine Soc Eur Spinal Deform Soc Eur Sect Cerv Spine Res Soc.* oct 2002;11 Suppl 2:S164-169.

16. S negas J, Vital JM, Pointillart V, Mangione P. Clinical evaluation of a lumbar interspinous dynamic stabilization device (the Wallis system) with a 13-year mean follow-up. *Neurosurg Rev.* juill 2009;32(3):335-41; discussion 341-342.
17. Gornet MF, Bailey AS, Eltahawy HA, Sasso RC, Cheng WK, Palmer S, et al. 250. Interspinous process distraction compared to nonoperative care for moderate lumbar degenerative disc disease: results of the FDA IDE clinical trial. *Spine J.* 1 sept 2019;19(9, Supplement):S122-3.
18. Sur YJ, Kong CG, Park JB. Survivorship analysis of 150 consecutive patients with DIAM™ implantation for surgery of lumbar spinal stenosis and disc herniation. *Eur Spine J Off Publ Eur Spine Soc Eur Spinal Deform Soc Eur Sect Cerv Spine Res Soc.* f vr 2011;20(2):280-8.
19. Daentzer D, Hurschler C, Seehaus F, Noll C, Schwarze M. Posterior dynamic stabilization in the lumbar spine - 24 months results of a prospective clinical and radiological study with an interspinous distraction device. *BMC Musculoskelet Disord.* 18 f vr 2016;17:90.
20. Jia YH, Sun PF. Preliminary evaluation of posterior dynamic lumbar stabilization in lumbar degenerative disease in Chinese patients. *Chin Med J (Engl).* janv 2012;125(2):253-6.
21. Deer T, Sayed D, Michels J, Josephson Y, Li S, Calodney AK. A Review of Lumbar Spinal Stenosis with Intermittent Neurogenic Claudication: Disease and Diagnosis. *Pain Med Malden Mass.* 1 d c 2019;20(Suppl 2):S32-44.
22. Sobottke R, Schl ter-Brust K, Kaulhausen T, R llinghoff M, Joswig B, St tzer H, et al. Interspinous implants (X Stop, Wallis, Diam) for the treatment of LSS: is there a correlation between radiological parameters and clinical outcome? *Eur Spine J Off Publ Eur Spine Soc Eur Spinal Deform Soc Eur Sect Cerv Spine Res Soc.* oct 2009;18(10):1494-503.
23. Zaina F, Tomkins-Lane C, Carragee E, Negrini S. Surgical Versus Nonsurgical Treatment for Lumbar Spinal Stenosis. *Spine.* 15 juill 2016;41(14):E857-68.
24. Zhao XW, Ma JX, Ma XL, Li F, He WW, Jiang X, et al. Interspinous process devices(IPD) alone versus decompression surgery for lumbar spinal stenosis(LSS): A systematic review and meta-analysis of randomized controlled trials. *Int J Surg Lond Engl.* mars 2017;39:57-64.
25. Deer TR, Sayed D, Malinowski MN, Rowe JJ, Jameson JB, Liang K, et al. A Review of Emerging Evidence for Utilization of a Percutaneous Interspinous Process Decompression Device to Treat Symptomatic Lumbar Adjacent-Segment Degeneration. *Pain Med Off J Am Acad Pain Med.* d c 2019;20(Suppl 2):S9-13.
26. Cairns K, Deer T, Sayed D, van Noort K, Liang K. Cost-effectiveness and Safety of Interspinous Process Decompression (Superion). *Pain Med Malden Mass.* 1 d c 2019;20(Suppl 2):S2-8.
27. Gazzeri R, Galarza M, Alfieri A. Controversies about interspinous process devices in the treatment of degenerative lumbar spine diseases: past, present, and future. *BioMed Res Int.* 2014;2014:975052.
28. Ge Y, Lu Y, Ma C, Lu B, Ma E, Zhang Y, et al. Effect of Different Interventions on Lumbar Spinal Stenosis: A Systematic Evaluation and Network Meta-Analysis. *World Neurosurg.* f vr 2025;194:123459.
29. Kaye AD, Edinoff AN, Temple SN, Kaye AJ, Chami AA, Shah RJ, et al. A Comprehensive Review of Novel Interventional Techniques for Chronic Pain: Spinal Stenosis and Degenerative Disc Disease—MILD Percutaneous Image Guided Lumbar Decompression, Vertiflex Interspinous Spacer, MinuteMan G3 Interspinous-Interlaminar Fusion. *Adv Ther.* 1 sept 2021;38(9):4628-45.
30. Aguwa OK, Park DK. Indirect spinal decompression for lumbar stenosis. *Semin Spine Surg.* 25. 2013;274-82.
31. Giraldo JP, Williams GP, Lee JJ, Potts EA, Uribe JS. An evidence-based review of the current surgical treatments for chronic low-back pain: rationale, indications, and novel therapies. *J*

- Neurosurg Spine. 1 févr 2025;1-12.
32. Senegas J, Vital JM, Pointillart V, Mangione P. Long -term actuarial survivorship analysis of an interspinous stabilization system. *Eur Spine J*. 2007;1279-89.
 33. Barrey CY, Le Huec JC, French Society for Spine Surgery. Chronic low back pain: Relevance of a new classification based on the injury pattern. *Orthop Traumatol Surg Res OTSR*. avr 2019;105(2):339-46.
 34. Marsh GDJ, Mahir S, Leyte A. A prospective randomised controlled trial to assess the efficacy of dynamic stabilisation of the lumbar spine with the Wallis ligament. *Eur Spine J Off Publ Eur Spine Soc Eur Spinal Deform Soc Eur Sect Cerv Spine Res Soc*. oct 2014;23(10):2156-60.
 35. Zini C, Bellini M, Masala S, Marcia S. Percutaneous Interspinous Spacer in Spinal-Canal-Stenosis Treatment: Pros and Cons. *Medicina (Mex)* [Internet]. 16 juill 2019 [cité 3 mars 2021];55(7). Disponible sur: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6681403/>
 36. Chen Z, Peng B, Li D, Pang X, Yang H. Minimum 5-year follow-up study on the effects of the Wallis dynamic stabilization system in the treatment of lumbar degenerative disease. *Chin Med J (Engl)*. 2014;127(20):3587-91.
 37. Gu H, Chang Y, Zeng S, Gao Y, Zhan S, Zheng Q. Efficacy of the Wallis interspinous implant for primary lumbar disc herniation : a prospective randomised controlled trial. *Acta Orthop Belg*. sept 2017;83(3):405-15.
 38. Gu H, Chang Y, Zeng S, Zheng X, Zhang R, Zhan S, et al. Wallis Interspinous Spacer for Treatment of Primary Lumbar Disc Herniation: Three-Year Results of a Randomized Controlled Trial. *World Neurosurg*. déc 2018;120:e1331-6.
 39. Korovessis P, Syrimpeis V, Tsekouras V, Vardakastanis K, Fennema P. PEEK versus Silicon Interspinous Spacer for Reduction of Supradjacent Segment Degeneration following Decompression and Short- Segment Instrumentation for Degenerative Lumbar Spinal Stenosis. *Adv Orthop* [Internet]. 8 août 2018 [cité 3 mars 2021];2018. Disponible sur: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6106718/>
 40. Zhou Z, Xiong W, Li L, Li F. Adjacent segmental degeneration following Wallis interspinous stabilization implantation. *Medicine (Baltimore)* [Internet]. 2 juin 2017 [cité 2 mars 2021];96(22). Disponible sur: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5459735/>
 41. Crawford RJ, Malone QJ, Price RI. Change of sagittal spinal alignment and its association with pain and function after lumbar surgery augmented with an interspinous implant. *Scoliosis Spinal Disord*. 2017;12:2.
 42. Jiang Y qi, Che W, Wang H ren, Li R yu, Li X lei, Dong J. Minimum 5 year follow-up of multi-segmental lumbar degenerative disease treated with discectomy and the Wallis interspinous device. *J Clin Neurosci Off J Neurosurg Soc Australas*. juill 2015;22(7):1144-9.
 43. Pan B, Zhang ZJ, Lu YS, Xu WG, Fu CD. Experience with the second-generation Wallis interspinous dynamic stabilization device implanted in degenerative lumbar disease: a case series of 50 patients. *Turk Neurosurg*. 2014;24(5):713-9.
 44. Roitberg B, Zileli M, Sharif S, Anania C, Fornari M, Costa F. Mobility-Preserving Surgery for Lumbar Spinal Stenosis: WFNS Spine Committee Recommendations. *World Neurosurg X*. juill 2020;7:100078.
 45. Wang H, Peng J, Zeng Q, Zhong Y, Xiao C, Ye Y, et al. Dynesys system vs posterior decompression and fusion for the treatment of lumbar degenerative diseases. *Medicine (Baltimore)*. 22 mai 2020;99(21):e19784.
 46. Machado GC, Ferreira PH, Yoo RI, Harris IA, Pinheiro MB, Koes BW, et al. Surgical options for

- lumbar spinal stenosis. Cochrane Database Syst Rev. 1 nov 2016;11:CD012421.
47. Yue ZJ, Liu RY, Lu Y, Dong LL, Li YQ, Lu EB. Middle-period curative effect of posterior lumbar intervertebral fusion (PLIF) and interspinous dynamic fixation (Wallis) for treatment of L4/5 degenerative disease and its influence on adjacent segment degeneration. *Eur Rev Med Pharmacol Sci.* déc 2015;19(23):4481-7.
 48. Cai Y, Luo J, Huang J, Lian C, Zhou H, Yao H, et al. Interspinous spacers versus posterior lumbar interbody fusion for degenerative lumbar spinal diseases: a meta-analysis of prospective studies. *Int Orthop.* juin 2016;40(6):1135-42.
 49. Li M, Yang H, Wang G. Interspinous process devices for the treatment of neurogenic intermittent claudication: a systematic review of randomized controlled trials. *Neurosurg Rev.* oct 2017;40(4):529-36.
 50. Poetscher AW, Gentil AF, Ferretti M, Lenza M. Interspinous process devices for treatment of degenerative lumbar spine stenosis: A systematic review and meta-analysis. *PloS One.* 2018;13(7):e0199623.
 51. Wei FL, Zhou CP, Liu R, Zhu KL, Du MR, Gao HR, et al. Management for lumbar spinal stenosis: A network meta-analysis and systematic review. *Int J Surg Lond Engl.* janv 2021;85:19-28.
 52. Ravindra VM, Ghogawala Z. Is There Still a Role for Interspinous Spacers in the Management of Neurogenic Claudication? *Neurosurg Clin N Am.* juill 2017;28(3):321-30.
 53. Filippiadis DK, Marcia S, Ryan A, Beall DP, Masala S, Deschamps F, et al. New Implant-Based Technologies in the Spine. *Cardiovasc Intervent Radiol.* oct 2018;41(10):1463-73.
 54. Ramesh A, Lyons F, Kelleher M. Aperius interspinous device for degenerative lumbar spinal stenosis: a review. *Neurosurg Rev.* avr 2016;39(2):197-205; discussion 205