



IMPLANT CARD

Instructions for completion





Symbols are explained on the back:









IMPLANT CARD

Instructions for completion

To be filled by the healthcare institution/provider:

- ①  Name of the patient or patient ID
- ②  Date of implantation
- ③  Name and address of the healthcare institution/provider
- ④  Manufacturer's information website

The Implant Card information corresponding to 5 to 10 are provided in form of one sticker with each device:

- ⑤  Device type in required language
- ⑥  Device name
- ⑦  UDI-DI Code (HRI)
- ⑧  Lot Number/ Batch Code
- ⑨  UDI code (AIDC format)
- ⑩  Name and address of the manufacturer of the implanted medical device

The device information sticker corresponding to each device used during surgery is available in the packaging and must be placed on each empty face of the implant card as shown below:



Per SPR 23.423.4, the Patient Leaflet inside the packaging of each device used during the operation contains all of the following particulars: (MDR Article 18 Text –Clause 1 – (b), (c), (d)). The Patient Leaflet of each device used should be supplied to the patient with the implant card.